










































 <p>VERMEIREN we care for you</p>				Ul. Łączna 1 55-100 Trzebnica Polska				
<h1>FOREST 3</h1> <p>Electric wheelchair with rear wheel drive. Excellent off-road properties.</p>				Delivery address				
Purchaser:						<u>ref.</u>		
<h2>Forest 3 standard version</h2>				Цена		4162.00 EUR		
КОНФИГУРАЦИЯ				ЦЕНИ в EUR с ДДС				
Motors								
	6 km/h				4162			
	10 km/h				4617			
	14 km/h				5021			
XD	PowerModule 120A - increases power - facilitates overcoming obstacles and driving on slopes				549			
SE 50	V-STAB TECHNOLOGY: automatic compensation of traction errors				660			
Seat width								
cm	40	45	50		-			
Seat depth								
cm	42	44	46	48	50	52	-	
Seat height (measured to the seat frame without the cushion)								
cm	41,5	43,5 ¹	46,5	48 ²	48,5 ³		-	
	¹ only in combination with the seat option SE43 (electric angle adjustment)							
	² only in combination with the seat option SE42 (electric elevator)							
	³ only in combination with the seat option SE43 (electric angle adjustment)							
Controller & Armrests								
Controller								
R	Right side				-			
L	Left side				-			
Tilting to the side								
SE53	Tiltable to the sides (to the outside and to the inside)				-			




















SE53Z	Tilted to the side and lowered under the armrest (inside and outside). Not compatible with SE64 - fixed & SE64 - swivel.		74	
SE53A	Tilts outwards and inwards with adjustable three-level locking. Not compatible with SE64 - fixed & SE64 - swivel.		40	
Controller options				
REM200	LINX controller (REM 200)		-	
REM400	LINX LCD controller with integrated LCD touch screen		550	
SE64 - fixed	Activation buttons - attached to the armrest. Not compatible with SE53Z & SE53A.		370	
SE64 - swivel	Activation buttons - attached to a tilt-up controller. Not compatible with SE53Z & SE53A.		385	
Special controllers				
SE68	Special levers for tetraplegics to switch functions on a LINX LCD joystick) - Only with REM 400 controller.		385	
SE02	Table with joystick (Requires LINX or LINX LCD driver to be ordered)		316	
SE04	<i>Chin control - Requires LINX LCD driver (REM400) to be ordered</i>		1410	
SE08	The controller is mounted behind the backrest for an accompanying person only		105	
SE09	Additional controller for accompanying person		610	
SE63	Foam ball on the controller		105	
SE82	<i>Control buttons located in various places (octopus 2). Requires a LINX LCD driver (REM400) to be ordered. And specify exactly where the buttons should be placed on the trolley.</i>		2020	
SE83 kulka	<i>Mini Joystick (minimal pressure force) - ball-shaped tip. Requires a LINX LCD driver (REM400) to be ordered.</i>		2690	
SE83 walek	<i>Mini Joystick (minimum pressure force) - a cylinder-shaped tip that is indented. Requires a LINX LCD driver (REM400) to be ordered.</i>		4230	
SE 55	<i>HD joystick - reinforced for people with high spasticity. Requires a LINX LCD driver (REM400) to be ordered.</i>		1670	
Armrest & Elbow Support				
	PU elbow support		-	
	Elbow support - Soft comfort		60	
	Support for the elbow - leather comfort		60	

B66	<i>PU gutter armrest (stabilizing the forearm), not compatible with B15 and B15A-tables, SE74A-mirror.</i>					
	Right		70			
	Left		70			
B66S	<i>Vacuum armrest, modified STABILO (stabilizing in front of the arm), not compatible with B15-table, SE74A-mirror.</i>					
	Right		175			
	Left		175			
Footrests						
B06	Removable double-plate footrests			-		
BA6	Removable double-plate footrests - shorter 7 cm			-		
BK6	removable single-plate footrests			60		
BL6	single-plate removable footrests - shorter 7 cm			60		
Manually adjustable footrests						
Right side						
BZ7	adjustable footrest to the level with calf supports			75		
BZ8	regulowana podnózek do poziomu z podparami pod łydki - poziom 7 cm			75		
Left side						
BZ7	adjustable footrest to the level with calf supports			75		
BZ8	regulowana podnózek do poziomu z podparami pod łydki - poziom 7 cm			75		
Electrical adjustment of the footrests						
Right side						
BZ7E	Adjustable footrest to the level with calf supports			550		
BZ8E	Adjustable footrest to the level with calf supports - shorter 7 cm			550		
Lewa strona						
BZ7E	Adjustable footrest to the level with calf supports			550		
BZ8E	Adjustable footrest to the level with calf supports - shorter 7 cm			550		
ADJUSTABLE HORIZONTAL FOOTRESTS WITH ONE FOOT AND CALF PLATE						
Bz10	Footrest with adjustable angle to horizontal, with one-piece calf support			150		
Seat (seat cushion fastened with Velcro straps to the seat plate L16)						
	Type of upholstery material					
	Nylon		-			
	PU Soft (grey)		-			
	AirFlow (black)		-			
	AirFlow (black) / cordura (black)		75			
	Cordura		60			
	Skaj (black)		-			










L16	plate with velcro stripes without cushion		-	
B23	4 cm flat pillow with velcro stripes		-	
B23P	Cushion made of foam in the cover		60	
L35	Contur cushion with velcro stripes		-	
VICAIR LIBERTY	VICAIR pneumatic cushion		120	
VICAIR PROFILE	VICAIR pneumatic cushion with pelvic stabilizer		170	
VICAIR VECTOR O2	8-chamber pneumatic cushion VICAIR VECTOR O2- 6 cm		340	
	8-chamber pneumatic cushion VICAIR VECTOR O2- 10 cm		340	
VICAIR MULTIFUNCTIONAL O2	5-chamber pneumatic cushion VICAIR MULTIFUNCTIONAL O2 - 6 cm		370	
VICAIR POMMEL O2	5-chamber pneumatic cushion VICAIR POMMEL O2- 10 cm		570	
VICAIR CENTRE RELIFE O2	5-chamber pneumatic cushion VICAIR CENTER RELIFE - 10 cm		395	
VICAIR ACTIVE O2	5-chamber pneumatic cushion VICAIR ACTIVE O2- 6 cm		345	
VICAIR ACTIVE O2	5-chamber pneumatic cushion VICAIR ACTIVE O2- 9 cm		345	
Backrest				
	Type of upholstery material			
	Nylon		-	
	PU Soft (grey)		-	
	AirFlow (black)		-	
	AirFlow (black) / cordura (black)		80	
	Cordura		60	
	Skaj (black)		-	
L34 Std	Cushion on a wooden base - Contour STD			
	Back height			
	58 cm		-	
	higher backrest 65 cm		55	
	higher backrest 70 cm			
higher backrest 75 cm				

L34 Tub	Tubular frame: cushion on a wooden base - Contour. With adjustable settings higher / lower and deeper / shallower.				
	Back height				
	58 cm			-	
	higher backrest 65 cm			55	
	higher backrest 70 cm				
	higher backrest 75 cm				
L34 Tub Ergo	Tubular frame: Ergo aluminum plate base. With adjustable settings higher / lower and deeper / shallower. (always a plate height of 53 cm).			65	
L44 Velcro	Tubular frame: elastic base with VELCRO straps			-	
Backrest cushions on L34 Tube Ergo board and L44 Velcro belts					
B23	Thin black nylon pillow - 4 cm - always in L44 Velcro.			-	
L34	Contur cushion			-	
Vicair Multifunctional	Vicair Multifunctional cushion			620	
BASE	STABILO backrest system - vacuum forming of the backrest.			375	
COMFORTABLE PLUS	STABILO backrest system - vacuum forming of the backrest.			430	
TARTA ° torso brace (only with L44 backrest and L34 Tube Ergo backrest) additional backrest to be mounted on the backrest frame.					
5S	5 vertebrae standard			1450	
5C	5 vertebrae contour			1450	
6S	6 vertebrae standard			1500	
6C	6 vertebrae contour			1500	
Seat & backrest adjustments					
SE30	Manual step adjustment of the seat angle			-	 SE34/SE36 SE39
SE34	Manual step adjustment of the backrest angle			-	
SE36	Manual, smooth adjustment of the backrest angle up to 30 ° - gas spring			240	
SE39	Automatic, smooth adjustment of the seat angle (20 °) and backrest (30 °) (2 actuators) - seat height adjustment (without a pillow) in the range of 41.5 cm or 46.5 cm.			1190	

SE40	Automatic smooth adjustment of the seat angle (20 ° - 1 actuator) - seat height adjustment (without a pillow) in the range of 41.5 cm or 46.5 cm.			605	
SE41	Automatic, smooth adjustment of the backrest angle (30 ° - 1 actuator) - seat height adjustment (without a pillow) in the range of 41.5 cm or 46.5 cm.			605	
SE42+SE40	The automatic lift raises the seat by 300 mm + automatic, smooth adjustment of the seat angle (45 °) (2 actuators) - seat height. (without pillow) 48 cm			1950	SE42+SE40/SE39
SE42+SE39	The automatic lift raises the seat by 300 mm + automatic, smooth adjustment of the seat angle (50 °) and backrest (30 °) (3 actuators) - seat height (without pillow) 48 cm			2200	
SE43	Automatic, smooth adjustment of the seat angle (50 °) (1 actuator) - seat height adjustment (without a pillow) in the range of 43.5 or 48.5 cm			620	
SE43+SE41	Automatic, smooth adjustment of the seat angle (50 °) and backrest (30 °) (2 actuators) - seat height adjustment. (without the pillow) in the range of 43.5 or 48.5 cm			1190	
Wheels & Tires					
Front wheels 9 '(230)					
	Black pneumatic tire, black rim, flat "line" tread			-	
	Pneumatic gray tire, black rim, flat "line" tread			-	
	Anti-puncture, black tire, black rim, flat "line" tread			105	
Front wheels 10 '(260)					
	Black pneumatic tire, black rim, thick block tread			-	
	Black pneumatic tire, black rim, flat "line" tread			-	
	Pneumatic gray tire, black rim, thick block tread			-	
	Anti-puncture black tire, black rim, thick block tread			130	
	Anti-puncture, black tire, black rim, flat "line" tread			130	
Rear driving wheel (3.00 x 8)					
	Black pneumatic tire, black rim, thick block tread			-	
	Black pneumatic tire, black rim, flat "line" tread			-	
	Pneumatic, gray tire, black rim, thick block tread			-	
	Anti-puncture, black tire, black rim, thick block tread			130	
Frame color					
			<i>matt</i>	<i>gloss</i>	
C4	Bright red				-
C75	Black matt			-	-
C77	Black gloss		-		-
C8	Yellow				-
C29	Silver				-
C59	Champagne				-
C53	Orange				-
					

C55	Inox			-	
C56	Blue metallic	-		-	
C76	White			-	
C79	Light blue			-	
C85	Light green			-	
C71	Navy blue shiny	-		-	
C72	Shiny red	-		-	
C73	Shiny lime green	-		-	
C78	Shiny bright blue	-		-	
C88	Yellow	-		-	
C89	Copper	-		-	
	Any color you choose			130	
Options					
Se 01	fenders			140	
B15	Plexiglass table (transparent) not compatible with the SE74 and SE74A mirror and the B66 and B66S armrests.			130	
B15A	<i>ABS table (black) not compatible with the SE74 and SE74A mirror and the B66 and B66S armrests. ATTENTION! Joystick must always be tilted to the side, Requires option SE53 to be ordered.</i>			75	
B22	Wedge			90	
B31	Crutch holder not compatible with SE42 & SE43 options			90	
B58	Safety belt with car buckle			-	
L04	Side pelotte stabilizing the torso (always with an L34 backrest on an upholstered wooden board)				
	Right big			75	
	Right small			75	
	Left big			75	
	Left small			75	
	Side pelotte stabilizing the torso (with standard L44 and L34 Tube Ergo backrest)				
	Right big			270	
Right small			270		

L04	Left big		270	
	Left small		270	
SE16	Side plate of the armrest (plexiglass) - transparent. Not compatible with SE37 femoral supports and B66 / B66S armrests.		-	
L55+B47	Orthopedic headrest adjustable in height and depth for standard L44, L34 Ergo upholstery.		140	
L55	Orthopedic headrest adjustable in height and depth in tandem with the L34 Std and L34 Tube backrest on an upholstered wooden board.		150	
L58+B47	Orthopedic 3D headrest, adjustable in height, depth, width and rotatable for standard L44, L34 Ergo upholstery.		180	
L58	Orthopedic 3D headrest, adjustable in height, depth, width and paired with the L34 Std and L34 Tube backrest on an upholstered wooden board.		185	
SE19	LED lights.		-	
SE58	"TAXI" handles.		-	
SE24	Bottle or cup holder.		60	
SE26	Push handle only with backrest option L34 Std.		76	
SE27	<i>Lockable trunk. Not compatible with L44 Velcro L34 Tube backrest and options, SE42 & SE43 & L55 & L58 & SE26.</i>		160	
SE28	Trunk		195	
SE37	Thigh pads (pair) - not compatible with the side plate of the SE16 armrest.		155	
SE74	<i>Rearview mirror not compatible with B66 / B66S armrests and B15 and B15A tables.</i>		55	
SE74A	<i>Tilting rear view mirror is not compatible with B66 / B66S armrests and B15 and B15A tables.</i>		110	
SE 06	USB port for charging the phone.		270	
SE 49	A bag for handy things under the seat, closed with a zipper.		60	
Belts, stabilizing vests				
U173 S	Torso stabilizing vest with a zipper		70	
U173 L			75	
U174 S	Torso stabilizing vest without a zipper		60	
U174 L			70	

U174 P-S	A vest stabilizing the torso with a zipper with an attachment to the crotch		60	
U174 P-L			70	
U180 S	Abdominal belt		45	
U180 L			50	
U180 P-S	Abdominal belt with a crotch opening		65	
U180 P-L			95	
U08 S	Abduction thigh belts		50	
U08 L			55	
U96	Wrist stabilization belts		35	
STABILIZING AND CORRECTION BELTS made of NEOFLEX ~ neoprene				
PELVIS				
U76 Neo	Two-point pelvic stabilizing belt		45	 
U75 Neo	A four-point pelvic stabilizing belt		45	
U01 Neo	Abduction belt		90	
TORSO				
U73 Neo	Torso stabilizing vest with a zipper		117	
U74 Neo	Torso stabilizing vest without a zipper		89	
U79 Neo	Four-point belts stabilizing the torso		76	
U80 Neo	A wide belt stabilizing the torso		62	
REMARKS!				

Reset the entire form.

The form will be sent to the sales department of Vermeiren Polska

VERMEIREN
we care for you