




# Normandie XXL

MAX  
users weight  
**350 kg**

## GERIATRIC XXL CHAIR

The Normandie is a geriatric chair with a reclinable backrest up to 45°. Equipped with a legrest that reclines automatically with the inclination of the backrest.

See page 88 for the options

-   
**C32**  
yellow
-   
**C33**  
light purple
-   
**C34**  
orange
-   
**C35**  
light green
-   
**C36**  
dark red
-   
**C37**  
night blue
-   
**C38**  
brown
-   
**C39**  
deep green
-   
**C40**  
blue fabric
-   
**C41**  
red fabric



Reclinable back 45°  
and legrest will recline  
automatically with the back

Armrests  
adjustable  
in height

# Altitude - Coraille XXL

MAX  
users weight  
**150 kg**

-   
dartex green
-   
dartex purple
-   
dartex praline
-   
dartex light purple

5 points safety  
harness

## GERIATRIC XXL CHAIR

Shell chair Altitude and Coraille XXL  
Seat width 550 mm or 600 mm for  
the overall width of 790 mm  
Maximum users weight: 150 kg

Only available with  
separated footrests



378 - Table on wheels



Model	890	740	0-220	510	1200	750	480	0-40°	780	38	350
<b>Normandie</b>	890	740	0-220	510	1200	750	480	0-40°	780	38	350
<b>Altitude / Coraille</b>	790	550	230	520.540.560	1240	1120	500	-4/36°	780	37	150
	790	600	230	520.540.560	1240	1120	500	-4/36°	780	37	150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

**XXL COMFORT**

