

o n e

Compact indoor scooter

Perfect for daily trips

Adapted to your style



o n e = 4

4 different colours included

Change your scooter's style whenever you want



VERMEIREN
we care for you

Compact and stylish

One is a perfect compact and stylish scooter for in- and outdoor use. Adapted to your style and your lifestyle! Lightweight and compact in its form, One can be used for short daily trips to restaurants, stores or in small spaces. The scooter fully disassembles for compact travelling or storage. 4 unique colours can be interchanged easily for a matching scooter to your style, at any time.

Colours

Your One scooter comes with 4 interchangeable colours included. Easy to change and adjustable to style. One = 4



Blue



Bright green



Orange



Soft Pink

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, **One** complies with the highest norms, laid down by certified test institutes in most European countries



Get on and off easily with the swivel seat



Fully demountable for compact storage

Technical information

■ Total width	490 mm
■ Total height	880 mm
■ Total length	1030 mm
■ Seat height	525-570 mm
■ Turning radius	1150 mm
■ Battery radius	24V - 15 Ah

■ Driving range	16 km
■ Motors	270 W
■ Speed	8 km/u
■ Weight	44,9 kg
■ Maximum user weight	136 kg
■ Extra informatie	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
België

www.vermeiren.be
info@vermeiren.be



For more information, contact your local dealer

VERMEIREN
we care for you